

EXPERIENCE STRENGTH & VERSATILITY


PRECISION SPINE
REFORM[®]
POCT
SYSTEM

Total freedom tulips
ALLOWS EASE OF ROD PLACEMENT

Smooth shank screw
PROTECTS NEURAL ELEMENTS

Complete set of connectors, dominoes & hooks
ENHANCES VERSATILITY

Low profile polyaxial screws
OFFERS 50° ANGULATION FOR INCREASED
INTRAOPERATIVE FLEXIBILITY



Midline screw fixation
PROVIDES STABLE FIXATION

Straight, lordosed & transitional rods
OPTIMIZES VERSATILITY

Rod-to-rod & tulip-to-tulip crosslinks
OFFERS 30° OF ROTATION AND
ANGULATION FOR SIMPLIFIED
PLACEMENT

Self-starting, aggressive screw tip
ALLOWS FOR IMMEDIATE BONE
ENGAGEMENT

PRECISION SPINE[®]
Discover the Difference



REFORM® POCT POLYAXIAL, SMOOTH SHANK & OCCIPITAL BONE SCREWS

Ø 3.5mm POLYAXIAL SCREWS

P/N	Description
55-PA-3508	3.5 x 8mm - POCT Screw
55-PA-3510	3.5 x 10mm - POCT Screw
55-PA-3512	3.5 x 12mm - POCT Screw
55-PA-3514	3.5 x 14mm - POCT Screw
55-PA-3516	3.5 x 16mm - POCT Screw
55-PA-3518	3.5 x 18mm - POCT Screw
55-PA-3520	3.5 x 20mm - POCT Screw
55-PA-3522	3.5 x 22mm - POCT Screw
55-PA-3524	3.5 x 24mm - POCT Screw
55-PA-3526	3.5 x 26mm - POCT Screw
55-PA-3528	3.5 x 28mm - POCT Screw
55-PA-3530	3.5 x 30mm - POCT Screw
55-PA-3535	3.5 x 35mm - POCT Screw
55-PA-3540	3.5 x 40mm - POCT Screw

Ø 4.0mm POLYAXIAL SCREWS

P/N	Description
55-PA-4008	4.0 x 8mm - POCT Screw
55-PA-4010	4.0 x 10mm - POCT Screw
55-PA-4012	4.0 x 12mm - POCT Screw
55-PA-4014	4.0 x 14mm - POCT Screw
55-PA-4016	4.0 x 16mm - POCT Screw
55-PA-4018	4.0 x 18mm - POCT Screw
55-PA-4020	4.0 x 20mm - POCT Screw
55-PA-4022	4.0 x 22mm - POCT Screw
55-PA-4024	4.0 x 24mm - POCT Screw
55-PA-4026	4.0 x 26mm - POCT Screw
55-PA-4028	4.0 x 28mm - POCT Screw
55-PA-4030	4.0 x 30mm - POCT Screw
55-PA-4035	4.0 x 35mm - POCT Screw
55-PA-4040	4.0 x 40mm - POCT Screw

Ø 4.5mm POLYAXIAL SCREWS

P/N	Description
55-PA-4520	4.5 x 20mm - POCT Screw
55-PA-4525	4.5 x 25mm - POCT Screw
55-PA-4530	4.5 x 30mm - POCT Screw
55-PA-4535	4.5 x 35mm - POCT Screw
55-PA-4540	4.5 x 40mm - POCT Screw
55-PA-4545	4.5 x 45mm - POCT Screw
55-PA-4550	4.5 x 50mm - POCT Screw

Ø 5.0mm POLYAXIAL SCREWS*

P/N	Description
55-PA-5020	5.0 x 20mm - POCT Screw
55-PA-5025	5.0 x 25mm - POCT Screw
55-PA-5030	5.0 x 30mm - POCT Screw
55-PA-5035	5.0 x 35mm - POCT Screw
55-PA-5040	5.0 x 40mm - POCT Screw
55-PA-5045	5.0 x 45mm - POCT Screw
55-PA-5050	5.0 x 50mm - POCT Screw

Ø 5.3mm POLYAXIAL SCREWS*

P/N	Description
55-PA-5320	5.3 x 20mm - POCT Screw
55-PA-5325	5.3 x 25mm - POCT Screw
55-PA-5330	5.3 x 30mm - POCT Screw
55-PA-5335	5.3 x 35mm - POCT Screw
55-PA-5340	5.3 x 40mm - POCT Screw
55-PA-5345	5.3 x 45mm - POCT Screw
55-PA-5350	5.3 x 50mm - POCT Screw

Ø 3.5mm SMOOTH SHANK SCREWS

P/N	Description
55-PB-3520	3.5 x 20mm - Smooth Shank
55-PB-3522	3.5 x 22mm - Smooth Shank
55-PB-3524	3.5 x 24mm - Smooth Shank
55-PB-3526	3.5 x 26mm - Smooth Shank
55-PB-3528	3.5 x 28mm - Smooth Shank
55-PB-3530	3.5 x 30mm - Smooth Shank
55-PB-3532	3.5 x 32mm - Smooth Shank
55-PB-3534	3.5 x 34mm - Smooth Shank

Ø 4.0mm SMOOTH SHANK SCREWS

P/N	Description
55-PB-4020	4.0 x 20mm - Smooth Shank
55-PB-4022	4.0 x 22mm - Smooth Shank
55-PB-4024	4.0 x 24mm - Smooth Shank
55-PB-4026	4.0 x 26mm - Smooth Shank
55-PB-4028	4.0 x 28mm - Smooth Shank
55-PB-4030	4.0 x 30mm - Smooth Shank
55-PB-4032	4.0 x 32mm - Smooth Shank
55-PB-4034	4.0 x 34mm - Smooth Shank

Ø 4.5mm OCCIPITAL SCREWS

P/N	Description
55-OS-4506	4.5 x 6mm - Plate Screw
55-OS-4508	4.5 x 8mm - Plate Screw
55-OS-4510	4.5 x 10mm - Plate Screw
55-OS-4512	4.5 x 12mm - Plate Screw
55-OS-4514	4.5 x 14mm - Plate Screw
55-OS-4516	4.5 x 16mm - Plate Screw

Ø 5.25mm OCCIPITAL SCREWS

P/N	Description
55-OS-5206	5.25 x 6mm - Plate Screw
55-OS-5208	5.25 x 8mm - Plate Screw
55-OS-5210	5.25 x 10mm - Plate Screw
55-OS-5212	5.25 x 12mm - Plate Screw
55-OS-5214	5.25 x 14mm - Plate Screw
55-OS-5216	5.25 x 16mm - Plate Screw

REFORM POCT SYSTEM TRAYS

P/N	Description
55-BK-0101	Standard Implants
55-BK-0102	Add-On Implants
55-BK-0103	Occhipital Implants
55-BK-0201	Standard Instruments
55-BK-0203	Add-On Instruments
55-BK-0202	Occhipital Instruments



55-BK-0201



55-BK-0203

REFORM POCT PLATES, RODS & CAP SCREW

OCCIPITAL PLATES

P/N	Description
55-OA-0101	Y-Shape Plate, 25-35mm
55-OA-0102	Y-Shape Plate, 32-42mm
55-OA-0103	Y-Shape Plate, 40-50mm
55-OA-0201	A-Shape Plate, 25-35mm
55-OA-0202	A-Shape Plate, 32-42mm
55-OA-0203	A-Shape Plate, 40-50mm

Ø 3.5mm STRAIGHT Ti RODS

P/N	Description
55-ST-5020	3.5 x 20mm - Rod
55-ST-5030	3.5 x 30mm - Rod
55-ST-5040	3.5 x 40mm - Rod
55-ST-5050	3.5 x 50mm - Rod
55-ST-5060	3.5 x 60mm - Rod
55-ST-5070	3.5 x 70mm - Rod
55-ST-5080	3.5 x 80mm - Rod
55-ST-5100	3.5 x 100mm - Rod
55-ST-5120	3.5 x 120mm - Rod
55-ST-5240	3.5 x 240mm - Rod**
55-TT-5502	250mm Ti Tran. Rod
55-TT-5508	550mm Ti Tran. Rod

LOCKING CAP

P/N	Description
55-LS-0100	Cap Screw (Lock-Screw)

OCCIPITAL Ti SHORTRODS

P/N	Description
55-OR-0500	3.5 x 90°
55-OR-0503	3.5 x 75°

OCCIPITAL CoCr SHORT RODS

P/N	Description
55-OR-0600	3.5 X 90°
55-OR-0603	3.5 X 75°

Ø 3.5mm LORDOTIC Ti RODS

P/N	Description
55-LT-5050	3.5 X 50mm*
55-LT-5060	3.5 x 60mm
55-LT-5070	3.5 X 50mm*
55-LT-5080	3.5 x 80mm
55-LT-5090	3.5 X 50mm*
55-LT-5100	3.5 x 100mm
55-LT-5050	3.5 X 110mm*
55-LT-5050	3.5 X 120mm*

STRAIGHT CoCr RODS

P/N	Description
55-SC-5080	3.5 x 80mm
55-SC-5120	3.5 x 120mm
55-SC-5240	3.5 x 240mm*
55-TC-5502	250mm CC Tran. Rod
55-TC-5508	550mm CC Tran. Rod

PRECISION SPINE®

www.precisionspineinc.com



Made in the USA

PRECISION SPINE
REFORM®
POCT
SYSTEM

*Available by Special Request
** Located in the Occipital Tray