

EXPERIENCE EFFICIENCY & VISIBILITY

PRECISION SPINE
REFORM[®]Ti
MODULAR TITANIUM PEDICLE SCREW

Audible click during tulip attachment
ENSURES THAT THE TULIP IS CONNECTED

Consistent attachment force
ENHANCES EASE OF ATTACHMENT

Triple lead thread
**INCREASES SPEED & EFFICIENCY OF SCREW DELIVERY
AND REDUCES FATIGUE**

30° of angulation in all planes
ALLOWS FOR INTRAOPERATIVE FLEXIBILITY



Standard and reduction tulips
PROVIDES TRUE VERSATILITY

Enlarged T25 drive feature
OFFERS STABILITY DURING INSERTION

Minimally tapered screw thread
**LOWERS SCREW INSERTION
TORQUE WHILE OPTIMIZING PULL OUT
STRENGTH**

Self starting aggressive screw tip
**OPTIMIZES INITIAL BONE
ENGAGEMENT**



PRECISION SPINE[®]

Discover the Difference



SYSTEM TRAYS

P/N	Description
39-BK-0541	Reform Ti Modular Implants
39-BK-0520	Reform Ti Modular Instruments Tray 1
39-BK-0202	Reform Ti Standard Instruments Tray 2



BONE SCREWS

Ø 4.5mm SCREWS

P/N	Description
39-SB-4525*	4.5 x 25mm - Ti Modular Shank
39-SB-4530*	4.5 x 30mm - Ti Modular Shank
39-SB-4535*	4.5 x 35mm - Ti Modular Shank
39-SB-4540*	4.5 x 40mm - Ti Modular Shank
39-SB-4545*	4.5 x 45mm - Ti Modular Shank

Ø 5.5mm SCREWS

P/N	Description
39-SB-5530	5.5 x 30mm - Ti Modular Shank
39-SB-5535	5.5 x 35mm - Ti Modular Shank
39-SB-5540	5.5 x 40mm - Ti Modular Shank
39-SB-5545	5.5 x 45mm - Ti Modular Shank
39-SB-5550	5.5 x 50mm - Ti Modular Shank

Ø 6.5mm SCREWS

P/N	Description
39-SB-6530	6.5 x 30mm - Ti Modular Shank
39-SB-6535	6.5 x 35mm - Ti Modular Shank
39-SB-6540	6.5 x 40mm - Ti Modular Shank
39-SB-6545	6.5 x 45mm - Ti Modular Shank
39-SB-6550	6.5 x 50mm - Ti Modular Shank
39-SB-6555	6.5 x 55mm - Ti Modular Shank

Ø 7.5mm SCREWS

P/N	Description
39-SB-7535	7.5 x 35mm - Ti Modular Shank
39-SB-7540	7.5 x 40mm - Ti Modular Shank
39-SB-7545	7.5 x 45mm - Ti Modular Shank
39-SB-7550	7.5 x 50mm - Ti Modular Shank
39-SB-7555	7.5 x 55mm - Ti Modular Shank

Ø 8.5mm SCREWS

P/N	Description
39-SB-8540*	8.5 x 40mm - Ti Modular Shank
39-SB-8545*	8.5 x 45mm - Ti Modular Shank
39-SB-8550*	8.5 x 50mm - Ti Modular Shank
39-SB-8555*	8.5 x 55mm - Ti Modular Shank
39-SB-8560*	8.5 x 60mm - Ti Modular Shank
39-SB-8570*	8.5 x 70mm - Ti Modular Shank
39-SB-8580*	8.5 x 80mm - Ti Modular Shank

Ø 9.5mm SCREWS

P/N	Description
39-SB-9560*	9.5 x 60mm - Ti Modular Shank
39-SB-9570*	9.5 x 70mm - Ti Modular Shank
39-SB-9580*	9.5 x 80mm - Ti Modular Shank

* Per Request Only



LOCKING CAP, CONNECTORS & RODS

TULIPS

P/N	Description
39-MT-0401	Reform Ti Modular Tulip
39-MT-0402*	Reform Ti Modular Reduction Tulip



LOCKING CAP

P/N	Description
39-LS-0100	Cap Screw (Lock-Screw)



Ø 5.5mm LORDOTIC Ti RODS

P/N	Description
39-LT-5035	Lordotic Ti Rod, x 35mm
39-LT-5040	Lordotic Ti Rod, x 40mm
39-LT-5045	Lordotic Ti Rod, x 45mm
39-LT-5050	Lordotic Ti Rod, x 50mm
39-LT-5055	Lordotic Ti Rod, x 55mm
39-LT-5060	Lordotic Ti Rod, x 60mm
39-LT-5065	Lordotic Ti Rod, x 65mm
39-LT-5070	Lordotic Ti Rod, x 70mm
39-LT-5075	Lordotic Ti Rod, x 75mm
39-LT-5080	Lordotic Ti Rod, x 80mm
39-LT-5090	Lordotic Ti Rod, x 90mm
39-LT-5100	Lordotic Ti Rod, x 100mm
39-LT-5110	Lordotic Ti Rod, x 110mm
39-LT-5120	Lordotic Ti Rod, x 120mm



CROSS CONNECTORS

P/N	Description
39-CC-0035	Cross Connector 35mm
39-CC-0040	Cross Connector 40mm
39-CC-0048	Cross Connector 48mm
39-CC-0066	Cross Connector 66mm



Ø 5.5mm STRAIGHT Ti RODS

P/N	Description
39-ST-5080	Straight Ti Rod, x 80mm
39-ST-5100	Straight Ti Rod, x 100mm
39-ST-5120	Straight Ti Rod, x 120mm
39-ST-5200	Straight Ti Rod, x 200mm
39-ST-5400	Straight Ti Rod, x 400mm

