


EXPERIENCE AN OPTIMAL ENVIRONMENT FOR FUSION



Plasma-sprayed Ti Surface
**PROVIDES ENHANCED STABILITY,
IMPROVED IMAGING, AND
INCREASED MIGRATION RESISTANCE**

Large graft window
**FOR MAXIMUM BIOLOGICS
TO FACILITATE FUSION**

Tantalum markers
OPTIMIZE FLUOROSCOPIC VISUALIZATION



Dual HA Coating & Plasma-Sprayed Ti
**PROVIDES OPTIMAL ENVIRONMENT
FOR OSSEOINTEGRATION**

Osteoconductive HA Coating
**FACILITATES FORMATION OF BONY
STRUCTURE & MAY AID IN MORE
STABLE FIXATION**

PEEK
**OFFERS A MODULUS OF ELASTICITY
SIMILAR TO BONE TO RESIST
SUBSIDENCE AND STRESS SHIELDING**

PRECISION SPINE®
Discover the Difference



SYSTEM TRAYS

P/N	Description
62-BK-0602	ShurFit ACIF 2C Implants
62-BK-01XX	ShurFit ACIF 2C Instruments



SHURFIT® ACIF 2C CAGES

14x12mm CAGES

P/N	Description
ACIF05-05P-2C	ACIF 2C Cage 14x12 - 5° x 5mm
ACIF05-06P-2C	ACIF 2C Cage 14x12 - 5° x 6mm
ACIF05-07P-2C	ACIF 2C Cage 14x12 - 5° x 7mm
ACIF05-08P-2C	ACIF 2C Cage 14x12 - 5° x 8mm
ACIF05-09P-2C	ACIF 2C Cage 14x12 - 5° x 9mm
ACIF05-10P-2C	ACIF 2C Cage 14x12 - 5° x 10mm
ACIF05-11P-2C	ACIF 2C Cage 14x12 - 5° x 11mm
ACIF05-12P-2C*	ACIF 2C Cage 14x12 - 5° x 12mm
ACIF10-05P-2C	ACIF 2C Cage 14x12 - 10° x 5mm
ACIF10-06P-2C	ACIF 2C Cage 14x12 - 10° x 6mm
ACIF10-07P-2C	ACIF 2C Cage 14x12 - 10° x 7mm
ACIF10-08P-2C	ACIF 2C Cage 14x12 - 10° x 8mm
ACIF10-09P-2C	ACIF 2C Cage 14x12 - 10° x 9mm
ACIF10-10P-2C	ACIF 2C Cage 14x12 - 10° x 10mm
ACIF10-11P-2C	ACIF 2C Cage 14x12 - 10° x 11mm
ACIF10-12P-2C*	ACIF 2C Cage 14x12 - 10° x 12mm

16x14mm CAGES

P/N	Description
45-05-05P-2C	ACIF 2C Cage 16x14 - 5° x 5mm
45-05-06P-2C	ACIF 2C Cage 16x14 - 5° x 6mm
45-05-07P-2C	ACIF 2C Cage 16x14 - 5° x 7mm
45-05-08P-2C	ACIF 2C Cage 16x14 - 5° x 8mm
45-05-09P-2C	ACIF 2C Cage 16x14 - 5° x 9mm
45-05-10P-2C	ACIF 2C Cage 16x14 - 5° x 10mm
45-05-11P-2C	ACIF 2C Cage 16x14 - 5° x 11mm
45-05-12P-2C*	ACIF 2C Cage 16x14 - 5° x 12mm
45-10-05P-2C	ACIF 2C Cage 16x14 - 10° x 5mm
45-10-06P-2C	ACIF 2C Cage 16x14 - 10° x 6mm
45-10-07P-2C	ACIF 2C Cage 16x14 - 10° x 7mm
45-10-08P-2C	ACIF 2C Cage 16x14 - 10° x 8mm
45-10-09P-2C	ACIF 2C Cage 16x14 - 10° x 9mm
45-10-10P-2C	ACIF 2C Cage 16x14 - 10° x 10mm
45-10-11P-2C	ACIF 2C Cage 16x14 - 10° x 11mm
45-10-12P-2C*	ACIF 2C Cage 16x14 - 10° x 12mm

*By Request Only



5°



10°

SIZERS, INSERTER & TAMP

14x12mm SIZERS

P/N	Description
ACIFS05	ACIF Sizer 14x12 - 5° & 10° x 5mm
ACIFS06	ACIF Sizer 14x12 - 5° & 10° x 6mm
ACIFS07	ACIF Sizer 14x12 - 5° & 10° x 7mm
ACIFS08	ACIF Sizer 14x12 - 5° & 10° x 8mm
ACIFS09	ACIF Sizer 14x12 - 5° & 10° x 9mm
ACIFS10	ACIF Sizer 14x12 - 5° & 10° x 10mm
ACIFS11	ACIF Sizer 14x12 - 5° & 10° x 11mm
ACIFS12	ACIF Sizer 14x12 - 5° & 10° x 12mm



INSERTER

P/N	Description
02-9001	Inserter



16x14mm SIZERS

P/N	Description
45-905	ACIF Sizer 16x14 - 5° & 10° x 5mm
45-906	ACIF Sizer 16x14 - 5° & 10° x 6mm
45-907	ACIF Sizer 16x14 - 5° & 10° x 7mm
45-908	ACIF Sizer 16x14 - 5° & 10° x 8mm
45-909	ACIF Sizer 16x14 - 5° & 10° x 9mm
45-910	ACIF Sizer 16x14 - 5° & 10° x 10mm
45-911	ACIF Sizer 16x14 - 5° & 10° x 11mm
45-912	ACIF Sizer 16x14 - 5° & 10° x 12mm

TAMP

P/N	Description
02-9002	Tamp



PRECISION SPINE®
www.precisionspineinc.com

