

EXPERIENCE STABLE FIXATION

Low profile titanium tulip
REDUCES RISK OF FACET IMPINGEMENT

30° of angulation in all planes
INCREASES INTRAOPERATIVE FLEXIBILITY

Hydroxyapatite coating
PROMOTES OSSEOINTEGRATION BETWEEN BONE AND IMPLANT WITH GREATER PULL OUT STRENGTH

Triple lead thread
INCREASES SPEED & EFFICIENCY OF SCREW DELIVERY AND REDUCES FATIGUE

PRECISION SPINE REFORM[®]TiHA


HA COATED TITANIUM PEDICLE SCREW SYSTEM

Enlarged T25 drive feature
PROVIDES STABILITY DURING INSERTION

Minimally tapered screw thread
LOWERS SCREW INSERTION TORQUE WHILE OPTIMIZING PULL OUT STRENGTH

Sterile packaging
PREVENTS IMPLANT CONTAMINATION

Self starting aggressive screw tip
OPTIMIZES INITIAL BONE ENGAGEMENT



PRECISION SPINE[®]

Discover the Difference



SYSTEM TRAYS

P/N	Description
39-BK-0731	Reform Ti HA Standard Implants
39-BK-0721	Reform Ti HA Standard Instrument Tray 1
39-BK-0202	Reform Ti HA Standard Instrument Tray 2



BONE SCREWS

Ø 4.5mm SCREWS

P/N	Description
39-HT-4525	4.5 x 25mm - Reform Ti HA Screw
39-HT-4530	4.5 x 30mm - Reform Ti HA Screw
39-HT-4535	4.5 x 35mm - Reform Ti HA Screw
39-HT-4540	4.5 x 40mm - Reform Ti HA Screw
39-HT-4545	4.5 x 45mm - Reform Ti HA Screw

Ø 5.5mm SCREWS

P/N	Description
39-HT-5530	5.5 x 30mm - Reform Ti HA Screw
39-HT-5535	5.5 x 35mm - Reform Ti HA Screw
39-HT-5540	5.5 x 40mm - Reform Ti HA Screw
39-HT-5545	5.5 x 45mm - Reform Ti HA Screw
39-HT-5550	5.5 x 50mm - Reform Ti HA Screw

Ø 6.5mm SCREWS

P/N	Description
39-HT-6530	6.5 x 30mm - Reform Ti HA Screw
39-HT-6535	6.5 x 35mm - Reform Ti HA Screw
39-HT-6540	6.5 x 40mm - Reform Ti HA Screw
39-HT-6545	6.5 x 45mm - Reform Ti HA Screw
39-HT-6550	6.5 x 50mm - Reform Ti HA Screw
39-HT-6555	6.5 x 55mm - Reform Ti HA Screw

Ø 7.5mm SCREWS

P/N	Description
39-HT-7535	7.5 x 35mm - Reform Ti HA Screw
39-HT-7540	7.5 x 40mm - Reform Ti HA Screw
39-HT-7545	7.5 x 45mm - Reform Ti HA Screw
39-HT-7550	7.5 x 50mm - Reform Ti HA Screw
39-HT-7555	7.5 x 55mm - Reform Ti HA Screw

Ø 8.5mm SCREWS

P/N	Description
39-HT-8540	8.5 x 40mm - Reform Ti HA Screw
39-HT-8545	8.5 x 45mm - Reform Ti HA Screw
39-HT-8550	8.5 x 50mm - Reform Ti HA Screw
39-HT-8555	8.5 x 55mm - Reform Ti HA Screw
39-HT-8560	8.5 x 60mm - Reform Ti HA Screw
39-HT-8570	8.5 x 70mm - Reform Ti HA Screw
39-HT-8580	8.5 x 80mm - Reform Ti HA Screw

PRECISION SPINE REFORM[®] Ti HA

HA COATED TITANIUM PEDICLE SCREW SYSTEM

LOCKING CAP, CROSS CONNECTORS & RODS

LOCKING CAP

P/N	Description
39-LS-0100	Cap Screw (Lock-Screw)



Ø 5.5mm STRAIGHT TI RODS

P/N	Description
39-ST-5080	Straight Ti Rod - 5.5mm x 80mm
39-ST-5100	Straight Ti Rod - 5.5mm x 100mm
39-ST-5120	Straight Ti Rod - 5.5mm x 120mm
39-ST-5200	Straight Ti Rod - 5.5mm x 200mm
39-ST-5400	Straight Ti Rod - 5.5mm x 400mm



CROSS CONNECTORS

P/N	Description
39-CC-0035	Cross Connector 35mm
39-CC-0040	Cross Connector 40mm
39-CC-0048	Cross Connector 48mm
39-CC-0066	Cross Connector 66mm



Ø 5.5mm LORDOTIC TI RODS

P/N	Description
39-LT-5035	Lordotic Ti Rod - 5.5mm x 35mm
39-LT-5040	Lordotic Ti Rod - 5.5mm x 40mm
39-LT-5045	Lordotic Ti Rod - 5.5mm x 45mm
39-LT-5050	Lordotic Ti Rod - 5.5mm x 50mm
39-LT-5055	Lordotic Ti Rod - 5.5mm x 55mm
39-LT-5060	Lordotic Ti Rod - 5.5mm x 60mm
39-LT-5065	Lordotic Ti Rod - 5.5mm x 65mm
39-LT-5070	Lordotic Ti Rod - 5.5mm x 70mm
39-LT-5075	Lordotic Ti Rod - 5.5mm x 75mm
39-LT-5080	Lordotic Ti Rod - 5.5mm x 80mm
39-LT-5090	Lordotic Ti Rod - 5.5mm x 90mm
39-LT-5100	Lordotic Ti Rod - 5.5mm x 100mm
39-LT-5110	Lordotic Ti Rod - 5.5mm x 110mm
39-LT-5120	Lordotic Ti Rod - 5.5mm x 120mm



PRECISION SPINE[®]

www.precisionspineinc.com

