

## EXPERIENCE OPTIMAL STABILITY




Large graft windows  
**OPTIMIZE BIOLOGICAL INGROWTH**

Central strut  
**PROVIDES A STABLE PLATFORM**

Tantalum markers  
**AID IN VISUALIZATION  
& POSITIONING**

Self-distracting bulleted leading edge  
**ALLOWS FOR EASY INSERTION**



Removal thread  
**ALLOWS REPOSITIONING OF IMPLANT**

Eight footprints with 0° & 7° Lordotic angles  
**ASSIST IN CORRECTING PATIENT'S  
SAGITTAL PROFILE**

Aggressive tooth pattern  
**MAINTAINS CAGE PLACEMENT**

**PRECISION SPINE®**  
*Discover the Difference*



## SYSTEM TRAYS

| P/N        | Description                     |
|------------|---------------------------------|
| 60-BK-5000 | ShurFit® LLIF 7° Cages & Trials |
| 60-BK-6000 | ShurFit® LLIF 0° Cages & Trials |
| 60-BK-2000 | Implantation Instruments        |
| 60-BK-7000 | Box Chisels & Trial Rasps       |



## ShurFit® LLIF 0° PEEK CAGES & TRIALS

### 0° CAGES (18W x L x H)

| P/N          | Description                   |
|--------------|-------------------------------|
| 23-PC1845-08 | LLIF 0° Cage - 18 x 45 x 08mm |
| 23-PC1845-10 | LLIF 0° Cage - 18 x 45 x 10mm |
| 23-PC1845-12 | LLIF 0° Cage - 18 x 45 x 12mm |
| 23-PC1845-14 | LLIF 0° Cage - 18 x 45 x 14mm |
| 23-PC1845-16 | LLIF 0° Cage - 18 x 45 x 16mm |
| 23-PC1850-08 | LLIF 0° Cage - 18 x 50 x 08mm |
| 23-PC1850-10 | LLIF 0° Cage - 18 x 50 x 10mm |
| 23-PC1850-12 | LLIF 0° Cage - 18 x 50 x 12mm |
| 23-PC1850-14 | LLIF 0° Cage - 18 x 50 x 14mm |
| 23-PC1850-16 | LLIF 0° Cage - 18 x 50 x 16mm |
| 23-PC1855-08 | LLIF 0° Cage - 18 x 55 x 08mm |
| 23-PC1855-10 | LLIF 0° Cage - 18 x 55 x 10mm |
| 23-PC1855-12 | LLIF 0° Cage - 18 x 55 x 12mm |
| 23-PC1855-14 | LLIF 0° Cage - 18 x 55 x 14mm |
| 23-PC1855-16 | LLIF 0° Cage - 18 x 55 x 16mm |
| 23-PC1860-08 | LLIF 0° Cage - 18 x 60 x 08mm |
| 23-PC1860-10 | LLIF 0° Cage - 18 x 60 x 10mm |
| 23-PC1860-12 | LLIF 0° Cage - 18 x 60 x 12mm |
| 23-PC1860-14 | LLIF 0° Cage - 18 x 60 x 14mm |
| 23-PC1860-16 | LLIF 0° Cage - 18 x 60 x 16mm |

### 0° CAGES (22W x L x H)

| P/N          | Description                   |
|--------------|-------------------------------|
| 23-PC2245-08 | LLIF 0° Cage - 22 x 45 x 08mm |
| 23-PC2245-10 | LLIF 0° Cage - 22 x 45 x 10mm |
| 23-PC2245-12 | LLIF 0° Cage - 22 x 45 x 12mm |
| 23-PC2245-14 | LLIF 0° Cage - 22 x 45 x 14mm |
| 23-PC2245-16 | LLIF 0° Cage - 22 x 45 x 16mm |
| 23-PC2250-08 | LLIF 0° Cage - 22 x 50 x 08mm |
| 23-PC2250-10 | LLIF 0° Cage - 22 x 50 x 10mm |
| 23-PC2250-12 | LLIF 0° Cage - 22 x 50 x 12mm |
| 23-PC2250-14 | LLIF 0° Cage - 22 x 50 x 14mm |
| 23-PC2250-16 | LLIF 0° Cage - 22 x 50 x 16mm |
| 23-PC2255-08 | LLIF 0° Cage - 22 x 55 x 08mm |
| 23-PC2255-10 | LLIF 0° Cage - 22 x 55 x 10mm |
| 23-PC2255-12 | LLIF 0° Cage - 22 x 55 x 12mm |
| 23-PC2255-14 | LLIF 0° Cage - 22 x 55 x 14mm |
| 23-PC2255-16 | LLIF 0° Cage - 22 x 55 x 16mm |
| 23-PC2260-08 | LLIF 0° Cage - 22 x 60 x 08mm |
| 23-PC2260-10 | LLIF 0° Cage - 22 x 60 x 10mm |
| 23-PC2260-12 | LLIF 0° Cage - 22 x 60 x 12mm |
| 23-PC2260-14 | LLIF 0° Cage - 22 x 60 x 14mm |
| 23-PC2260-16 | LLIF 0° Cage - 22 x 60 x 16mm |

### 0° TRIALS (18W x H)

| P/N        | Description               |
|------------|---------------------------|
| 60-TS-0808 | LLIF 0° Trial - 18 x 08mm |
| 60-TS-0810 | LLIF 0° Trial - 18 x 10mm |
| 60-TS-0812 | LLIF 0° Trial - 18 x 12mm |
| 60-TS-0814 | LLIF 0° Trial - 18 x 14mm |
| 60-TS-0816 | LLIF 0° Trial - 18 x 16mm |

### 0° TRIALS (22W x H)

| P/N        | Description               |
|------------|---------------------------|
| 60-TS-0208 | LLIF 0° Trial - 22 x 08mm |
| 60-TS-0210 | LLIF 0° Trial - 22 x 10mm |
| 60-TS-0212 | LLIF 0° Trial - 22 x 12mm |
| 60-TS-0214 | LLIF 0° Trial - 22 x 14mm |
| 60-TS-0216 | LLIF 0° Trial - 22 x 16mm |

## ShurFit® LLIF 7° PEEK CAGES & TRIALS

### 7° CAGES (18W x L x H)

| P/N          | Description                   |
|--------------|-------------------------------|
| 23-EC1845-08 | LLIF 7° Cage - 18 x 45 x 08mm |
| 23-EC1845-10 | LLIF 7° Cage - 18 x 45 x 10mm |
| 23-EC1845-12 | LLIF 7° Cage - 18 x 45 x 12mm |
| 23-EC1845-14 | LLIF 7° Cage - 18 x 45 x 14mm |
| 23-EC1845-16 | LLIF 7° Cage - 18 x 45 x 16mm |
| 23-EC1850-08 | LLIF 7° Cage - 18 x 50 x 08mm |
| 23-EC1850-10 | LLIF 7° Cage - 18 x 50 x 10mm |
| 23-EC1850-12 | LLIF 7° Cage - 18 x 50 x 12mm |
| 23-EC1850-14 | LLIF 7° Cage - 18 x 50 x 14mm |
| 23-EC1850-16 | LLIF 7° Cage - 18 x 50 x 16mm |
| 23-EC1855-08 | LLIF 7° Cage - 18 x 55 x 08mm |
| 23-EC1855-10 | LLIF 7° Cage - 18 x 55 x 10mm |
| 23-EC1855-12 | LLIF 7° Cage - 18 x 55 x 12mm |
| 23-EC1855-14 | LLIF 7° Cage - 18 x 55 x 14mm |
| 23-EC1855-16 | LLIF 7° Cage - 18 x 55 x 16mm |
| 23-EC1860-08 | LLIF 7° Cage - 18 x 60 x 08mm |
| 23-EC1860-10 | LLIF 7° Cage - 18 x 60 x 10mm |
| 23-EC1860-12 | LLIF 7° Cage - 18 x 60 x 12mm |
| 23-EC1860-14 | LLIF 7° Cage - 18 x 60 x 14mm |
| 23-EC1860-16 | LLIF 7° Cage - 18 x 60 x 16mm |

### 7° CAGES (22W x L x H)

| P/N          | Description                   |
|--------------|-------------------------------|
| 23-EC2245-08 | LLIF 7° Cage - 22 x 45 x 08mm |
| 23-EC2245-10 | LLIF 7° Cage - 22 x 45 x 10mm |
| 23-EC2245-12 | LLIF 7° Cage - 22 x 45 x 12mm |
| 23-EC2245-14 | LLIF 7° Cage - 22 x 45 x 14mm |
| 23-EC2245-16 | LLIF 7° Cage - 22 x 45 x 16mm |
| 23-EC2250-08 | LLIF 7° Cage - 22 x 50 x 08mm |
| 23-EC2250-10 | LLIF 7° Cage - 22 x 50 x 10mm |
| 23-EC2250-12 | LLIF 7° Cage - 22 x 50 x 12mm |
| 23-EC2250-14 | LLIF 7° Cage - 22 x 50 x 14mm |
| 23-EC2250-16 | LLIF 7° Cage - 22 x 50 x 16mm |
| 23-EC2255-08 | LLIF 7° Cage - 22 x 55 x 08mm |
| 23-EC2255-10 | LLIF 7° Cage - 22 x 55 x 10mm |
| 23-EC2255-12 | LLIF 7° Cage - 22 x 55 x 12mm |
| 23-EC2255-14 | LLIF 7° Cage - 22 x 55 x 14mm |
| 23-EC2255-16 | LLIF 7° Cage - 22 x 55 x 16mm |
| 23-EC2260-08 | LLIF 7° Cage - 22 x 60 x 08mm |
| 23-EC2260-10 | LLIF 7° Cage - 22 x 60 x 10mm |
| 23-EC2260-12 | LLIF 7° Cage - 22 x 60 x 12mm |
| 23-EC2260-14 | LLIF 7° Cage - 22 x 60 x 14mm |
| 23-EC2260-16 | LLIF 7° Cage - 22 x 60 x 16mm |

### 7° TRIALS (18W x H)

| P/N        | Description               |
|------------|---------------------------|
| 60-TS-7808 | LLIF 7° Trial - 18 x 08mm |
| 60-TS-7810 | LLIF 7° Trial - 18 x 10mm |
| 60-TS-7812 | LLIF 7° Trial - 18 x 12mm |
| 60-TS-7814 | LLIF 7° Trial - 18 x 14mm |
| 60-TS-7816 | LLIF 7° Trial - 18 x 16mm |

### 7° TRIALS (22W x H)

| P/N        | Description               |
|------------|---------------------------|
| 60-TS-7208 | LLIF 7° Trial - 22 x 08mm |
| 60-TS-7210 | LLIF 7° Trial - 22 x 10mm |
| 60-TS-7212 | LLIF 7° Trial - 22 x 12mm |
| 60-TS-7214 | LLIF 7° Trial - 22 x 14mm |
| 60-TS-7216 | LLIF 7° Trial - 22 x 16mm |



PRECISION SPINE®

www.precisionspineinc.com

