

PRECISION SPINE

REFORM[®] Ti

MODULAR TITANIUM PEDICLE SCREW

VERSATILITY & VISIBILITY

The Reform[®] Ti Modular Pedicle Screw System provides surgeons with increased flexibility, versatility and visibility to meet the varying requirements of degenerative and trauma procedures. The titanium tulips are available in standard and reduction options, with an audible click to confirm attachment. The titanium modular screws feature a minimally tapered, triple lead thread as well as a self-starting aggressive screw tip and enlarged T25 drive feature — all of which combines to deliver more immediate, secure bone engagement and maximum control during insertion. The Reform Ti Modular intraoperatively assembled construct has a pull off strength as strong as a Reform preassembled screw, providing optimal security.



TITANIUM TULIP



REDUCTION TULIPS



TRIPLE LEAD THREAD



REFORM[®] Ti MODULAR PEDICLE SCREW

Self Starting Aggressive Screw Tip optimizes initial bone engagement



Standard & Reduction Tulips provide true versatility



Triple Lead Thread increases speed & efficiency of screw delivery and reduces fatigue



Audible Click During Tulip Attachment ensures that the tulip is connected



Minimally Tapered Screw Thread lowers screw insertion torque while optimizing pull out strength



Consistent Attachment Force enhances ease of attachment

30° Of Angulation In All Planes allows for intraoperative flexibility

Enlarged T25 Drive Feature offers stability during insertion



Reform Ti Modular Sizing

4.5mm Diameter	5.5mm Diameter	6.5mm Diameter	7.5mm Diameter	8.5mm Diameter	9.5mm Diameter
25mm	30mm	30mm	35mm	40mm	60mm
30mm	35mm	35mm	40mm	45mm	70mm
35mm	40mm	40mm	45mm	50mm	80mm
40mm	45mm	45mm	50mm	55mm	
45mm	50mm	50mm	55mm	60mm	
		55mm		70mm	
				80mm	

*Per Request Only

Reform Rod Sizing – ø5.5mm

Ti Straight	Ti Curved	CoCr Straight	CoCr Curved
80mm	35-80mm (5mm)	80mm	35-80mm (5mm)
100mm	90-120mm (10mm)	100mm	90-120mm (10mm)
120mm		120mm	
200mm		200mm	
300mm		300mm	
400mm		400mm	
500mm		500mm	

Precision Spine, Inc.
 2050 Executive Drive, Pearl, MS 39208
 Customer Service: 1.888.241.4773
 Phone: 601.420.4244
 Toll Free: 877.780.4370
 Fax: 601.420.5501
www.precisionspineinc.com



CE 2797