

PRECISION SPINE
SHURFIT[®]
INTERBODY CAGES


PATIENT MATCHING OPTIONS

The ShurFit family of interbody cages includes designs for T-PLIF, PLIF and TLIF procedures. The ShurFit designs offer a large contact area to optimize vertebral body support and minimize possible risk of subsidence. The large graft window allows for maximum biological coverage area for fusion. Simple instrumentation offers ease of implantation for surgeons, while the unique tooth pattern geometry may minimize the potential for expulsion. The PEEK ShurFit Interbody cages are available in multiple sizes to accommodate varying anatomies and effect optimal patient matching.

ShurFit[®] T-PLIF



ShurFit[®] PLIF




ShurFit[®] TLIF



AGGRESSIVE TOOTH PATTERN



BULLETED NOSE




LARGE BONE GRAFT WINDOW



SHURFIT® INTERBODY CAGES

ShurFit® PLIF Interbody Cages




- Large graft window
- Convex & Lordotic Series (Available in 5°, 10° & 15°)
- Bulleted Nose
- 22mm & 25mm Lengths
- 9mm Width
- 8 to 16mm Heights (1mm increments)

ShurFit® TLIF Interbody Cages



- Large graft window
- 24mm & 30mm Lengths
- 9mm Width
- 6 to 16mm Heights (2mm increments)

ShurFit® T-PLIF Interbody Cages



- Large graft window
- Straight and Convex Options
- Bulleted Nose
- 28.5mm Lengths
- 11mm Width
- 7 to 16mm Heights (1mm increments)

Three options from the extended line of ShurFit Interbody Cages.

Precision Spine, Inc.

2050 Executive Drive, Pearl, MS 39208

Customer Service: 1.888.241.4773

Phone: 601.420.4244

Toll Free: 877.780.4370

Fax: 601.420.5501

www.precisionspineinc.com

Caution: Federal (USA) law restricts these devices to sale by, or on the order of, a physician.

ShurFit® and Precision Spine® are trademarks of Precision Spine, Inc.

Copyright 2020 Precision Spine, Inc. All rights reserved. P/N LBL-FB-002 Rev E 03/2020

