


PRECISION SPINE
SLIMPLICITY[®] HP
ANTERIOR CERVICAL PLATING SYSTEM

SIMPLE AND VERSATILE

The Slimplicity[®] HP ACP System offers simplified versatility to accommodate diverse patient anatomies and pathologies. Low profile constrained, semi-constrained, or hybrid constructs are achieved via fixed and variable angle screws.

- **Dual thread screws** maximize cortical and cancellous bone interface
- **One-step locking mechanism** provides tactile and visual confirmation
- **Alignment Markings** and **Large Graft Windows** for increased intra-operative visualization and positioning



VARIABLE ANGULATION

23°



CORTICAL CANCELLOUS SCREWS



ONE STEP LOCKING




PRECISION SPINE[®]

Discover the Difference




SLIMPLICITY® HP ANTERIOR CERVICAL PLATING SYSTEM

VARIABLE




FIXED




BONE SCREWS

Variable



Fixed




Lengths



PLATE SIZING*

- **Width** 16.5mm
- **Waist** 13mm
- **Thickness** 2mm



*Measured Center Hole to Center Hole


PRECISION SPINE SHURFIT[®] ACIF 2C HA-CPTI COATED INTERBODY CAGES



- Ti Plasma Coated Surface
- Hydroxyapatite (HA) Coating
- Increased Expulsion Resistance
- Tri-Bond Lock Technology
- Enhanced Imaging Properties

SIZING

- 14mm Wide x 12mm Deep
- 16mm Wide x 14mm Deep
- Heights - 5-12mm
- 5° & 10° Lordotic Options



Precision Spine, Inc.
2050 Executive Drive, Pearl, MS 39208
Customer Service: 1.888.241.4773
Phone: 601.420.4244
Toll Free: 877.780.4370
Fax: 601.420.5501
www.precisionspineinc.com

