





PRODUCT SUMMARY




PRECISION SPINE
VaultTM
ALIF SYSTEM

Footprints Available



Set Screws



- Attaches the VAULT™ Plate to the VAULT Peek Cage
- Used with the VAULT 2.5mm driver
- T8 Hexalobe



Screws

Self-Drilling

- Diameter 4.0mm
- Titanium
- 10° Medial/lateral angulation
- 40° Cephalad (head) to caudal (tail) angulation
- Screw thread washout: 1mm
- Variable angle
- 2.5mm hex



4.0mm Diameter


Lengths

25mm
30mm
35mm




Install Guide Block

- Attaches to the VAULT Plate to act as guide for the Awl and Bone Screws
- Designed to absorb and distribute force during implant insertion
- Available in 11mm, 13mm, 15mm, 17mm and 19mm heights
- Color coded with the VAULT Locking Cap




Locking Cap

- Reduces risk of disengagement/back out of Bone Screws
- Available in 11mm, 13mm, 15mm, 17mm and 19mm heights
- Color coded with the Install Guide Blocks and matching the VAULT plate/cage heights




VAULT™ - Top Tray (Tray Number 21-1040-2-CA)



#	Part No.	Description	Qty
1	10-C3908-11	VAULT ALIF Peek Cage 39mm x 11mm 8 deg lordosis	3
2	10-C3908-13	VAULT ALIF Peek Cage 39mm x 13mm 8 deg lordosis	3
3	10-C3908-15	VAULT ALIF Peek Cage 39mm x 15mm 8 deg lordosis	3
4	10-C3908-17	VAULT ALIF Peek Cage 39mm x 17mm 8 deg lordosis	3
5	10-C3908-19	VAULT ALIF Peek Cage 39mm x 19mm 8 deg lordosis	3
6	10-IB30-11	VS ALIF Cage Plate Install Block 30 x 11	1
7	10-IB30-13	VS ALIF Cage Plate Install Block 30 x 13	1
8	10-IB30-15	VS ALIF Cage Plate Install Block 30 x 15	1
9	10-IB30-17	VS ALIF Cage Plate Install Block 30 x 17	8
10	10-IB30-19	VS ALIF Cage Plate Install Block 30 x 19	1
11	10-C3915-11	VAULT ALIF Peek Cage 39mm x 11mm 15 deg lordosis	3
12	10-C3915-13	VAULT ALIF Peek Cage 39mm x 13mm 15 deg lordosis	3
13	10-C3915-15	VAULT ALIF Peek Cage 39mm x 15mm 15 deg lordosis	3

#	Part No.	Description	Qty
14	10-C3915-17	VAULT ALIF Peek Cage 39mm x 17mm 15 deg lordosis	3
15	10-C3915-19	VAULT ALIF Peek Cage 39mm x 19mm 15 deg lordosis	3
16	10-P3908-11P	VAULT ALIF Peek Plate 39mm x 11mm 8 deg lordosis	3
17	10-P3908-13P	VAULT ALIF Peek Plate 39mm x 13mm 8 deg lordosis	3
18	10-P3908-15P	VAULT ALIF Peek Plate 39mm x 15mm 8 deg lordosis	3
19	10-P3908-17P	VAULT ALIF Peek Plate 39mm x 17mm 8 deg lordosis	3
20	10-P3908-19P	VAULT ALIF Peek Plate 39mm x 19mm 8 deg lordosis	3
21	10-P3915-11P	VAULT ALIF Peek Plate 39mm x 11mm 15 deg lordosis	3
22	10-P3915-13P	VAULT ALIF Peek Plate 39mm x 13mm 15 deg lordosis	3
23	10-P3915-15P	VAULT ALIF Peek Plate 39mm x 15mm 15 deg lordosis	3
24	10-P3915-17P	VAULT ALIF Peek Plate 39mm x 17mm 15 deg lordosis	3
25	10-P3915-19P	VAULT ALIF Peek Plate 39mm x 19mm 15 deg lordosis	3


VAULT™ - Bottom Tray (Tray Number 21-1040-2-CA)



#	Part No.	Description	Qty
1	10-P3208-11P	VAULT ALIF Peek Plate 32mm x 11mm 8 deg lordosis	3
2	10-P3208-13P	VAULT ALIF Peek Plate 32mm x 13mm 8 deg lordosis	3
3	10-P3208-15P	VAULT ALIF Peek Plate 32mm x 15mm 8 deg lordosis	3
4	10-P3208-17P	VAULT ALIF Peek Plate 32mm x 17mm 8 deg lordosis	3
5	10-P3208-19P	VAULT ALIF Peek Plate 32mm x 19mm 8 deg lordosis	3
6	10-1019	Slap Hammer	1
7	10-P3215-11P	VAULT ALIF Peek Plate 32mm x 11mm 15 deg lordosis	3
8	10-P3215-13P	VAULT ALIF Peek Plate 32mm x 13mm 15 deg lordosis	3
9	10-P3215-15P	VAULT ALIF Peek Plate 32mm x 15mm 15 deg lordosis	3
10	10-P3215-17P	VAULT ALIF Peek Plate 32mm x 17mm 15 deg lordosis	3
11	10-P3215-19P	VAULT ALIF Peek Plate 32mm x 19mm 15 deg lordosis	3

#	Part No.	Description	Qty
12	10-C3208-11	VAULT ALIF Peek Cage 32mm x 11mm 8 deg lordosis	3
13	10-C3208-13	VAULT ALIF Peek Cage 32mm x 13mm 8 deg lordosis	3
14	10-C3208-15	VAULT ALIF Peek Cage 32mm x 15mm 8 deg lordosis	3
15	10-C3208-17	VAULT ALIF Peek Cage 32mm x 17mm 8 deg lordosis	3
16	10-C3208-19	VAULT ALIF Peek Cage 32mm x 19mm 8 deg lordosis	3
17	10-C3215-11	VAULT ALIF Peek Cage 32mm x 11mm 15 deg lordosis	3
18	10-C3215-13	VAULT ALIF Peek Cage 32mm x 13mm 15 deg lordosis	3
19	10-C3215-15	VAULT ALIF Peek Cage 32mm x 15mm 15 deg lordosis	3
20	10-C3215-17	VAULT ALIF Peek Cage 32mm x 17mm 15 deg lordosis	3
21	10-C3215-19	VAULT ALIF Peek Cage 32mm x 19mm 15 deg lordosis	3


VAULT™ - Top Tray (Tray Number 21-1041-2-CA)



#	Part No.	Description	Qty
1	10-S3208-11	VS Sizer 8 degree 32mm x 11mm	1
2	10-S3208-13	VS Sizer 8 degree 32mm x 13mm	1
3	10-S3208-15	VS Sizer 8 degree 32mm x 15mm	1
4	10-S3208-17	VS Sizer 8 degree 32mm x 17mm	1
5	10-S3208-19	VS Sizer 8 degree 32mm x 19mm	1
6	10-S3908-11	VS Sizer 8 degree 39mm x 11mm	1
7	10-S3908-13	VS Sizer 8 degree 39mm x 13mm	1
8	10-S3908-15	VS Sizer 8 degree 39mm x 15mm	1
9	10-S3908-17	VS Sizer 8 degree 39mm x 17mm	1
10	10-S3908-19	VS Sizer 8 degree 39mm x 19mm	1
11	10-S3215-11	VS Sizer 15 degree 32mm x 11mm	1
12	10-S3215-13	VS Sizer 15 degree 32mm x 13mm	1
13	10-S3215-15	VS Sizer 15 degree 32mm x 15mm	1
14	10-S3215-17	VS Sizer 15 degree 32mm x 17mm	1
15	10-S3215-19	VS Sizer 15 degree 32mm x 19mm	1

#	Part No.	Description	Qty
16	10-S3915-11	VS Sizer 15 degree 39mm x 11mm	1
17	10-S3915-13	VS Sizer 15 degree 39mm x 13mm	1
18	10-S3915-15	VS Sizer 15 degree 39mm x 15mm	1
19	10-S3915-17	VS Sizer 15 degree 39mm x 17mm	1
20	10-S3915-19	VS Sizer 15 degree 39mm x 19mm	1
21	10-1001	VAULT Set Screw	8
22	10-1002-11	VAULT Locking Cap 11mm	4
23	10-1002-13	VAULT Locking Cap 13mm	4
24	10-1002-15	VAULT Locking Cap 15mm	4
25	10-1002-17	VAULT Locking Cap 17mm	4
26	10-1002-19	VAULT Locking Cap 19mm	4
27	10-4025	VAULT ALIF Self Drilling Variable Screw	8
28	10-4030	VAULT ALIF Self Drilling Variable Screw	8
29	10-4035	VAULT ALIF Self Drilling Variable Screw	8

VAULT™ - Bottom Tray (Tray Number 21-1041-2-CA)



#	Part No.	Description	Qty
1	10-1003-2	Universal Driver 2.5mm Rev 0 (PSSRS Handle)	2
2	10-1004-2	ALIF Cage Plate Inserter (Used with Install Guide Block)	2
3	10-1006-2	ALIF Cage Plate Sizer Inserter	2
4	10-1007-2	Fixed Awl 30 degree	1

#	Part No.	Description	Qty
5	10-1021	T-8 Hexalobe	2
6	10-1017	VAULT Removal Instrument	1
7	10-1014	Torque Driver 2.5mm Hex (Handle 10-1013)	1
8	10-1008-2	ALIF Cage Plate Universal Awl	1



Precision Spine, Inc.
2050 Executive Drive, Pearl, MS 39208
Customer Service: 1-888.241.4773
Fax: 601.420.5501

www.precisionspineinc.com

Please refer to package insert (LBL-FU-007) for complete system description, indications and warnings.
Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician.
Precision Spine™ and VAULT™ are trademarks of Precision Spine, Inc.
Copyright © 2015 Precision Spine, Inc. all rights reserved. LBL-FB-037 Rev B.